

# TABLETOP CENTRIFUGE CEN11-05



# TABLETOP CENTRIFUGE CEN11-05

These centrifuges are widely used in the fields of biochemistry, pharmacy, clinical medicine and other research fields for variety of separation needs. The programmable microprocessor centrifuge allows modification of the parameter during operation. In addition, accommodates various type of tubes or plates and size of the sample.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

## CEN11-05 TABLETOP CENTRIFUGE

Gentle braking at low speeds

LCD display

Short Spin key for quick centrifugation

Lid is opened automatically after the operation is finished, which avoids overheating of the sample and save processing time

Safety lid-lock



## SPECIFICATIONS

|                                       |                            |
|---------------------------------------|----------------------------|
| Model                                 | CEN11-05                   |
| Maximum Capacity (No of tubes x Vol.) | 6x15ml                     |
| Maximum Speed RPM                     | 100-5000                   |
| Maximum RCF                           | 3074xg                     |
| Timer Range                           | 15 S~99 min                |
| Acceleration/ Deceleration            | 20-80s↑ 20-80s↓            |
| Power                                 | 40 W                       |
| Power Supply                          | 110V 50/60Hz, 230V-50/60Hz |

## OPTIONAL ACCESSORIES

| Accessory Code | Name  | Description                 | Capacity |
|----------------|-------|-----------------------------|----------|
| LS55158        | Rotor | D12x75 glass tube           | 6x15ml   |
| LS55159        | Rotor | Conical tube with screw cap | 8x15 ml  |
| LS55160        | Rotor | Conical tube with screw cap | 12x5 ml  |



***Centrifugen***

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)