

LED OR SHADOWLESS LAMP OPL12-16



LED OR SHADOWLESS LAMP OPL12-16

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

OPL12-16 LED OR SHADOWLESS LAMP

Imported LED cold light source is adopted , with the ultrathin optical lens

Digital means are adopted for stepless regulation of LED brightness


Light field size can be adjusted by the handle

The removable handle sheath can be disinfected

The service life of LED shadowless lamp is long (50,000h)

The imported switch power supply is adopted to control the voltage, making the work voltage permanently stable

The no welding titanium alloy arm, Imported spring arm make sure the lamp head adjusted easily and freely



SPECIFICATIONS

Model	OPL12-16
Illumination	120000 Lux
Depth of light beam	≥500 mm
The value of the total irradiance	<6mW / m ²
Total irradiance	≤600W/m ²
Color rendering Index (Ra)	100≥Ra≥85
Color Temperature (K)	4000±500 K
Power supply	220±10% 50 Hz DC24 V
Power of bulb	300 W
Input power	350 VA

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com