

INCUBATOR SHAKER SHA2A-0300



INCUBATOR SHAKER SHA2A-0300

In Incubator Shaker, the constant agitation with uniform amplitude provides the optimum oxygen transfer to the cultures, leading to better and consistent yield. A variety of platforms and different laboratory accessories permits configuration to fulfill your demands.

Used in Cell Cultures, Cell Aeration, Microbiology, Increasing Solubility Rates, Metabolism Studies, Bacterial Cultures, Bacteriology.

Also known as Laboratory Incubator Shaker, Incubator Shaker.

SHA2A-0300 INCUBATOR SHAKER

The electronic control of closed cycle adjustable heating system

Sainless steel interior and Electrostatic sprayed exterior

Call to restore function, without power supply

With strong influence of intermittent rapid heating system of.200 type and 100B type with air pressure regulating valve

Flexible and convenient, wide variety of detachable accessories



SPECIFICATIONS

| | |
|------------------------|-----------------------------------|
| Model | SHA2A-0300 |
| Speed Range | 20-300 rpm |
| Standard Capacity | 500mlx4, 250mlx4, 100mlx5, 50mlx5 |
| Temperature | 4°C~60°C |
| Temperature Accuracy | ±0.1°C |
| Temperature Uniformity | ±1°C |
| Shaking Speed Accuracy | ±1 rpm |
| Timing Range | Optional |
| Motion | Orbital |
| Amplitude | 26 mm |
| Tray | 1 |
| Platform Dimension | 450x400 mm |
| Inner Dimension | 508x458x360 mm |
| Overall Dimension | 800x700x530 mm |
| Display | LCD |
| Weight | 100 kg |
| Power | 580 W |
| Power Supply | 220V±10% |

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com