

COAGULATION ANALYZER



COAGULATION ANALYZER

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods
Used in Clinical diagnostic, Laboratory, Pharmaceutical.
Also known as coagulation analyzer.

ANA13 AUTO COAGULATION ANALYZER



SPECIFICATIONS

| Model | ANA13-140 | ANA13-160 |
|------------------------|--------------------|---|
| Throughput | 140T/H | 160T/H |
| Analysis Method | Coagulation Method | Coagulation method, Immunospectrophotometry |
| Analysis Items | PT, APTT, TT, FIB | PT, APTT, TT, FIB, IDD, FDP |
| Sample Positions | 5 | 6 |
| Tilt reagent Positions | 6 | 11 refrigerated reagent positions |



ANA13-140



ANA13-160

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com