

STEREO ZOOM MICROSCOPE



STEREO ZOOM MICROSCOPE

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

MSC41 STEREO ZOOM MICROSCOPE


Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Convenient operation and special design for viewing incubating cell tissue.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



Product Image Coming Soon

SPECIFICATIONS

| Model | MSC41-045M | MSC41-045MT | MSC41-020ST |
|--|------------|-------------|-------------|
| Viewing Head | | | |
| Binocular view head inclined at 60° | Yes | No | - |
| Trinocular view head inclined at 45° | No | Yes | - |
| Digital view head inclined at 45° | No | Yes | - |
| Eyepiece | | | |
| High point wide field adjustable eyepiece WF10x22 mm | Yes | | - |
| MSC41-020ST | - | - | - |
| MSC41-040ST | - | - | - |
| MSC41-060ST | - | - | - |
| Zoom Objective | | | |
| 0.67 X—4.5 X | Yes | | - |
| Zoom Ratio | | | |
| 6.7:1 | Yes | | - |
| Auxiliary Objective | | | |
| 0.5X /177 mm, 1.5X /47 mm, 2X/26 mm | Yes | | - |
| Working Distance | | | |
| 115 mm | Yes | | - |
| MSC41-020ST | - | - | - |
| MSC41-040ST | - | - | - |
| MSC41-060ST | - | - | - |
| Focusing Bracket | | | |
| Focusing frame and adjustment knob | Yes | | - |
| Interpupillary Distance | | | |
| 55-75 mm / 52-75 mm | Yes | | - |

| | | | |
|---------------------------------|----------------|---|-------------------------|
| Up Illumination | | | |
| LED lamp | Yes | | - |
| Down Illumination | | | |
| LED lamp | Yes | | - |
| Power Supply | | | |
| AC110/220 V \pm 10%, 50/60 Hz | Yes | | - |
| MSC41-020ST | - | - | - |
| MSC41-040ST | - | - | - |
| MSC41-060ST | - | - | - |
| Package | | | |
| Cartoon with foam | Yes | | - |
| Package Size | 500x260x370 mm | | 300x200x150 mm |
| Gross Weight | 8 kg | | 3 kg |
| Illumination | - | - | Natural light Binocular |
| View Head And Body | - | - | Straight Head |
| Objective | - | - | 2X |

| Model | MSC41-040ST | MSC41-060ST |
|--|-------------|-------------|
| Viewing Head | | |
| Binocular view head inclined at 60° | - | - |
| Trinocular view head inclined at 45° | - | - |
| Digital view head inclined at 45° | - | - |
| Eyepiece | | |
| High point wide field adjustable eyepiece WF10x22 mm | - | - |
| MSC41-020ST | - | - |
| MSC41-040ST | - | - |
| MSC41-060ST | - | - |
| Zoom Objective | | |
| 0.67 X—4.5 X | - | - |
| Zoom Ratio | | |
| 6.7:1 | - | - |
| Auxiliary Objective | | |
| 0.5X /177 mm, 1.5X /47 mm, 2X/26 mm | - | - |
| Working Distance | | |
| 115 mm | - | - |
| MSC41-020ST | - | - |
| MSC41-040ST | - | - |
| MSC41-060ST | - | - |
| Focusing Bracket | | |
| Focusing frame and adjustment knob | - | - |
| Interpupillary Distance | | |
| 55-75 mm / 52-75 mm | - | - |
| Up Illumination | | |
| LED lamp | - | - |
| Down Illumination | | |
| LED lamp | - | - |
| Power Supply | | |
| AC110/220 V \pm 10%, 50/60 Hz | - | - |

| | | |
|--------------------|------------------------------------|--|
| MSC41-020ST | - | - |
| MSC41-040ST | - | - |
| MSC41-060ST | - | - |
| Package | | |
| Cartoon with foam | - | - |
| Package Size | 360x170x290 mm | 400x250x360 mm |
| Gross Weight | 6 kg | 8 kg |
| Illumination | Tranilluminator / inclined light | Natural light |
| View Head And Body | Inclined viewing Up-rotatable body | Inclined viewing Body 360°rotatable |
| Objective | 2X&4X, 90°rotated | 2X&4X, 360°rotated |



Product Image Coming Soon

MSC41-045M



MSC41-045MT



MSC41-020ST



MSC41-040ST



MSC41-060ST

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com